


# Digital Abbe Refractometer



## Technical Data

## WYA-2S

Measuring Range	1.3000-1.7000 nD 0 ~ 95% Brix
Accuracy	nD: $\leq \pm 0.0002$ Brix: $\leq \pm 0.1\%$
Temp. Resolution	0.1°C
Correction Range Brix v/s Temp.	15~45°C
Usage Temperature Range	35°C (room temperature)
Interface	RS232
Lamp	6.5v, 0.3A
Power Source	AC 220V / 50Hz

## RAVI SCIENTIFIC INDUSTRIES

14-UA Jawahar Nagar, Delhi - 110007, Delhi, India

Phone: 011-23857757 / 23855156

Mob/Whatsapp: 9811319003, 9910475789

[www.raviscientific.in](http://www.raviscientific.in)

[sales@raviscientific.in](mailto:sales@raviscientific.in)

