


Silicone Tubing (Food Grade)



- Translucent natural color for visual contact with the flow
- Odorless, tasteless, and inert

- Non-yellowing, soft & flexible
- Resists continuous temperature of 200°C



ID = Inner Diameter
OD = Outer Diameter
WT = Wall Thickness

Silicone may be low pressure steam sterilized in-line or autoclaved at up to 121°C in a normal autoclaving cycle. However, if exposed to repeated steam sterilization or long-term high temperature or pressure, Silicone will eventually relax and become gummy.

Model	Inner Dia	Outer Dia	Wall Thickness
RST-46	4	x 6	1
RST-69	6	x 9	1.5
RST-610	6	x 10	2
RST-810	8	x 10	1
RST-812	8	x 12	2
RST-1014	10	x 14	2
RST-1217	12	x 17	2.5

RAVI SCIENTIFIC INDUSTRIES

14-UA Jawahar Nagar, Delhi - 110007, Delhi, India

Phone: 011-23857757 / 23855156

Mob/Whatsapp: 9811319003, 9910475789

www.raviscientific.in

sales@raviscientific.in

