

Hollow Cathode Lamp



- High spectral response
- Stable light output

- Low noise characteristics
- Standard and encoded lamps



Technical Data

Diameter	1.5" (37mm) & 2" (50mm)
Functions	Standard and encoded lamps Long life due to higher gas volume Clear signals and low distortion due to special cathode technology Can be used for almost any AA instrument Low noise characteristics
Models Nos	Suffix 'HCL' with the atomic number Example: Atomic number of Zinc is 30 So Zinc Hollow Cathode Lamp's model number becomes HCL-30.

RAVI SCIENTIFIC INDUSTRIES

14-UA Jawahar Nagar, Delhi - 110007, Delhi, India

Phone: 011-23857757 / 23855156

Mob/Whatsapp: 9811319003, 9910475789

www.raviscientific.in

sales@raviscientific.in

