

# Digital Viscometer



- Large, backlit LCD display
- With durable carrying case

- High accuracy microprocessor
- Includes 4 spindle rotors



## Technical Data

## **A-220B**

Measuring Range	20 ~ 2,000,000 (mPa.s)
No. of Speeds	0.3, 0.6, 1.5, 3, 6, 12, 30, 60 RPM (8 speeds)
Spindles	#1, #2, #3, #4 is standard configuration
Accuracy	±1.0% Full Scale
Repeatability	±0.5% Full Scale
Display	LCD
Temp. Reading	Yes
Functions	With auto range function Liner Calibration Sound Alarm at under 20% Torque Time Function High quality Stepping Motor
Power Source	100V-240V

## **RAVI SCIENTIFIC INDUSTRIES**

14-UA Jawahar Nagar, Delhi - 110007, Delhi, India

Phone: 011-23857757 / 23855156

Mob/Whatsapp: 9811319003, 9910475789

[www.raviscientific.in](http://www.raviscientific.in)

[sales@raviscientific.in](mailto:sales@raviscientific.in)

